

[Document Name] Specification

[Title of the Invention] NOVEL PSYCHROTROPHIC BACTERIUM AND
DNA PROBE FOR DETECTING THE BACTERIUM

[Claim]

[Claim 1] A 16S rDNA which has the base sequence of SEQ ID NO:1.

[Claim 2] An oligonucleotide probe which comprises part of the base
sequence of SEQ ID NO:1.

[Claim 3] The oligonucleotide probe according to claim 2 wherein said
part of the base sequence of SEQ ID NO:1 comprises the base sequence of
SEQ ID NO:2.

[Claim 4] The oligonucleotide probe according to claim 2 or 3 for
detecting or identifying a bacterium selected from the group consisting of
Psychrobacter pacificus, *Psychrobacter glacincola* and analogs thereof.

[Claim 5] A method for detecting or identifying a bacterium selected
from the group consisting of *Psychrobacter pacificus*, *Psychrobacter
glacincola* and analogs thereof, using an oligonucleotide probe comprising
part of the base sequence of SEQ ID NO:1.

[Claim 6] *Psychrobacter pacificus*, which is aerobic, gram-negative,
nonmotile, colorless, non-sporulating and oxidase-positive.

[Detailed Description of the Invention]

[0001]

[Field of the Invention]

The present invention relates to a technique for monitoring the
circulation and upwelling of deep-sea water using a deep-sea
microorganism as an indicator, and particularly to a technique for
species-specifically detecting a bacterium selected from a group consisting
of *Psychrobacter pacificensis*, *Psychrobacter glacincola* and analogs
thereof.

[0002]

[Prior Art]

It is thought that waters in the deep-sea near the Japan islands are supplied by a great ocean current which originally starts from the high-density sea water in the depth of the sea near Greenland and via the Antarctic region where other high-density sea water merges with the current, flows to northern parts of the Pacific Ocean including the Japan Trench. Such deep-sea water contains plenty of nutritive salts and exhibits a high productivity of organisms in its upwelling area. Accordingly, the industrial applicability of such deep-sea waters is now explored.

[0003]

Further, deep-sea fish, which had not been utilized to date, are now beginning to be used as food or feed.

Moreover, it is suggested in connection with problems of global or local environmental pollution, that CO₂, radioactive waste or industrial waste resulting from human activities be disposed into the deep-sea zone near Japan.

[0004]

However, because there is little knowledge concerning deep-sea waters or zones, it is difficult to evaluate what effect deep-sea water has on the activities of organisms inhabiting in epipelagic zones, or what effect the disposal of, for example, CO₂, radioactive waste or industrial waste into the deep-sea would have on the activities of organisms in the deep-sea. Further, no indicator organism has been reported which may provide useful information on the global ocean current of deep-layer sea water.

[0005]

[Problems to be Solved by the Invention]

The object of the present invention is to provide a technique for evaluating the biological safety of artificial use of deep-sea water or deep-sea

zones, and particularly a technique for species-specifically detecting a microorganism naturally inhabiting in the deep sea or an analog thereof, based on the characteristics of its genetic information.

[0006]

[Means for Solving the Problem]

The present inventors developed an oligonucleotide probe which enables species-specific detection of a novel psychrotrophic bacterium species isolated from the deep-sea water of the Japan Trench, at a molecular or cell level based on base sequence information of 16S rRNA derived from the microorganism, thereby completing the present invention.

[0007]

In summary, the present invention provides a 16S rDNA which has the base sequence of SEQ ID NO:1.

The present invention also provides an oligonucleotide probe which comprises part of the base sequence of SEQ ID NO:1. The oligonucleotide probe may be either an RNA or DNA probe. One example of said part of the base sequence of SEQ ID NO:1 is the base sequence of SEQ ID NO:2. Such a probe can be used to detect or identify a bacterium selected from the group consisting of *Psychrobacter pacificus*, *Psychrobacter glacincola* and analogs thereof.

[0008]

The present invention also provides a method for detecting or identifying a bacterium selected from the group consisting of *Psychrobacter pacificus*, *Psychrobacter glacincola* and analogs thereof, using an oligonucleotide probe comprising part of the base sequence of SEQ ID NO:1.

Further, the present invention provides *Psychrobacter pacificus* which is aerobic, gram-negative, nonmotile, colorless, non-sporulating and oxidase-positive.

[0009]

[Embodiments of the Invention]

The newly-discovered microorganism species according to the present invention, *Psychrobacter pacificus*, is a heterotrophic microorganism which predominantly appears under cold culture conditions of 1 atm, at 4°C, that was isolated from the seawater of the Japan Trench in the offing of Hachijo island, Japan at a depth of 6,000 meters. Six strains of *Psychrobacter pacificus*, NIBH P2J2, NIBH P2J3, NIBH P2J13, NIBH P2K6, NIBH P2K17 and NIBH P2K18 have been isolated by the present inventors. Among these strains, phylogenetic characteristics of NIBH P2J3, NIBH P2K6 and NIBH P2K18 are shown in Tables 1-3 below. Tables 1-3 also include other strains isolated simultaneously with, and analogs of, the above-described strains, together with their phylogenetic characteristics.

[0010]

Table 1
Motility and Extracellular Organ of Psychrotrophic Bacteria Isolated from
Surface and Deep Seawaters of the Japan Trench

Strain	Motility test ¹ (microscopic)	Motility test ² (agar plate)	Extracellular organ ³	Phylogenetic location
Surface seawater				
P1H8	-	-	Flagella*	<i>Halomonadaceae</i>
P1H13	-	-	None	<i>Halomonadaceae</i>
P1H14	-	-	None	<i>Halomonadaceae</i>
P1H22	+	+	Flagella	<i>Halomonadaceae</i>
P1H25	+	+	Flagella	<i>Halomonadaceae</i>
Deep seawater				
P2J2	-/w	-	Fimbriae	<i>Moraxellaceae</i>
P2J3	-/w	-	Fimbriae	<i>Moraxellaceae</i>
P2J13	-/w	-	Fimbriae	<i>Moraxellaceae</i>
P2K6	-/w	-	Fimbriae	<i>Moraxellaceae</i>
P2K17	-/w	-	Fimbriae	<i>Moraxellaceae</i>
P2K18	-/w	-	Fimbriae	<i>Moraxellaceae</i>

[0011]

1: By optical microscopy using Nomarski optics. 2: On a semisolid agar medium with nutrient gradient. 3: By electron microscopy:- (negative); Flagella* (frequent adhesion of flagella was observed); w (weak twitch).

[0012]

Table 2
Characteristics of Phenotypes and GC Contents of
***Psychrobacter pacificensis* Strains and Their Analogs**

Characteristics ^a	<i>Psychrobacter pacificensis</i>				<i>Psychrobacter</i>	<i>Psychrobacter</i>	<i>Psychrobacter</i>	<i>Psychrobacter</i>	<i>Psychrobacter</i>
	NIBH strain no.				<i>immobilis</i> ^b	<i>urativorans</i> ^b	<i>frigidicola</i> ^b	<i>phenylpyruvicus</i>	<i>glaciacola</i>
	P2J3	P2K6	P2K18	Summary	(Phenon 1)	(Phenon 2)	(Phenon 3)	ACAM 535 ^b	ACAM 483 ^c
Urease activity	+	+	+	+	v+	v+	-	+	v-
Phenylalanine deaminase	-	-	-	-	+	-	+	+	-
Tryptophan deaminase	-	-	-	-	v-	-	+	-	-
Nitrate reduction	-	-	-	-	v-	v-	-	-	v+
Growth in NaCl (%):									
0	w	-	-	-	+	+	+	?	+
1	+	-	-	-v	+	+	+	+	+
3	+	+	+	+	+	+	+	+	+
5	+	+	+	+	+	+	+	+	+
8	-	-	-	-	+	+	+	+	+
Growth at (°C):									
30	+	+	+	+	+	-	-	+	-
35	+	w	+	+	v+	-	-	+	-
40	-	-	-	-	-	-	-	-	-
Acid production from:									
Glucose	+	w	+	+	-	-	-	-	-
Xylose	+	w	+	+	+	-	-	-	-
Arabinose	+	+	+	+	+	-	-	-	-
*Others ^d	-	-	-	-	-	-	-	-	-
PNPG test ^d	-	-	-	-	-	-	-	-	-
Use as sole carbon and energy sources:									
Acetate	+	-	-	-v	+	v+	+	+	+
L-alanine	+	-	-	-v	+	-	-	+	v+
p-hydroxy-benzoate	-	-	-	-	-	v-	-	-	?
3-hydroxy-butyrate	+	-	+	+v	+	+	-	+	v+
Citrate	-	-	-	-	v-	-	-	+	v+
Gluconate ^d	-	-	-	-	v-	-	-	-	-
L-histidine	+	+	+	+	+	-	-	-	v+
Lactate	+	-	+	+v	+(DL)	v+(DL)	-(DL)	+(DL)	v+(DL)
DL-malate ^d	+	+	+	+	+(L)	v+(L)	+(L)	+(L)	-(L)
Malonate	-	+	-	-v	-	-	-	-	-
Propionate	-	-	-	-	v+	-	-	+	+
L-serine	-	-	-	-	-	v-	-	-	-
Suberate	+	-	-	-v	v-	-	+	-	?
n-valerate	-	-	-	-	+	v+	+	+	+
**Others	-	-	-	-	-	-	-	-	-
DNA G+C content (mol %)	44	44	45	44-45	44-47	44-46	41-42	43	43-44

[0013]

a) All of the species and strains were proved to be positive for oxidase, catalase, culture at 4-15°C, resistance against 6.5% NaCl, as well as for the use of L-proline as sole carbon and energy sources. b) Data from Bowman et al., (1996) Int. J. Syst. Bacteriol. 46:841-848. c) Data from Bowman et al., (1997) System. Appl. Microbiol. 20:209-215. d) Determined by API 20 NE test. Compound availability was estimated by using API ID 32 GN test. PNPG is a test for β -galactosidase using para-nitrophenyl-(β)-D-galactopyranoside.

[0014]

*Others: glucose fermentation, indole production, hydrolysis of esculin, hydrolysis of gelatin, and arginine dihydrolase. ^(d)

**Others: N-acetyl-D-glucosamine, m-hydroxy-benzoate, glycogen, phenyl acetate. The following carbon and energy sources were not utilized by any species or strains: N-acetylglucosamine, adipic acid, L-arabinose, capric acid, L-fucose, 2-keto-gluconate, 5-keto-gluconate, (D)glucose, (myo) inositol, itaconic acid, maltose, D-mannitol, D-melibiose, (L) rhamnose, D-ribose, (D) salicin, D-sorbitol, sucrose. The type of substrate used by Bowman et al. (1996) is indicated in parenthesis. The type of optical isomer is indicated in parentheses in the Table.

Frequencies of positive strains in the columns of *Psychrobacter pacificensis* in Table 2 above: + = 100%; +v = 67%; -v = 33%; and - = 0%. Frequencies of positive strains in the columns of other *Psychrobacter* species in Table 2 above: + = 100-90%; v+ = 89-11%; and v- = 10-0%. w: weak response.

Strain NIBH P2K6 was defined as the reference strain of *Psychrobacter pacificensis*.

[0015]

Table 3

Fatty Acid Composition and Major Quinone Type of *Psychrobacter*
pacificensis

Composition	<i>Psychrobacter pacificensis</i>			Average content	<i>Psychrobacter immobilis</i> ATCC 43116
	NIBH strain no.				
	P2J3	P2K6*	P2K18		
Fatty acid:					
10:0	1.3	Tr	1.2	0.8 (0.7)	0.9
11:0	0.1	Tr	0.2	0.1 (0.1)	Tr
12:0	2.2	0.8	2.3	1.8 (0.8)	Tr
14:0	0.7	0.6	0.5	0.6 (0.1)	0.3
14:1 (X1)	0.1	0.2	Tr	0.1 (0.1)	0.1
15:0	0.4	Tr	0.4	0.3 (0.2)	0.2
16:0	7.3	8.7	6.5	7.5 (1.1)	4.3
16:1 (w7c)	9.7	15.8	6.6	10.7 (4.7)	3.8
16:1 (X2)	0.4	0.4	0.3	0.4 (0.1)	0.4
17:0	2.7	5.6	4.8	4.4 (1.5)	4.2
i17:0	0.6	Tr	0.4	0.3 (0.3)	Tr
17:1 (X3)	5.1	1.5	5.1	3.9 (2.1)	6.8
18:0	9.4	5.6	13.1	9.4 (3.8)	8.0
18:1 (w7c)	1.2	0.8	0.7	0.9 (0.3)	0.8
18:1 (w9c)	50.1	52.8	50.9	51.3 (1.4)	63.1
18:2	3.4	2.4	1.8	2.5 (0.8)	3.5
19:0	0.4	0.4	0.6	0.5 (0.1)	0.9
20:0	Tr	Tr	Tr	Tr	0.1
Total	95.1	95.6	95.4	95.5	97.4
Total unsaturated	70.0	73.9	65.4	69.8	75.0
Hydroxy fatty acid:					
3-OH 12:0	4.1	3.6	3.9	3.9 (0.3)	2.2
3-OH 14:0	0.8	0.8	0.7	0.8 (0.1)	0.4
Total	4.9	4.4	4.6	4.6 (0.3)	2.6
Major quinone type	Q-8	Q-8	Q-8	Q-8	Q-8

[0016]

X1-3: Accurate location of double bond has not been determined. Tr: trivial (<0.1%).

Numbers in parentheses indicate standard deviations (n=3). *: Reference strain.

Psychrobacter pacificus is an aerobic, gram-negative, nonmotile, colorless, non-sporulating, and oxidase-positive coccobacillus of 1.0-1.5 μ m long \times about 1 μ m wide. *Psychrobacter pacificensis* strains produce a number of fimbriae as extracellular organs but not flagella. Off-white colored, round, convex colonies with entire margins form on an agar plate containing polypeptone and yeast extract. No fluorescent color is formed. For optimum growth, seawater, or about 3% NaCl, may be required though most of the strains will not grow in the presence of 0%, 8% or higher NaCl. It takes 1-2 weeks for those strains to reach the stationary phase at 4°C, though these strains exhibit growth rates at 4°C comparable to their growth yields at 20°C. Optimum growth may be obtained at about 25°C, and the critical growth temperature is 38°C. Acids may be aerobically produced from glucose, xylose and arabinose. These strains are urease-activity-positive but phenylalanine deaminase- and tryptamine deaminase- negative. This species is negative in biochemical tests for glucose fermentation, indole production, hydrolysis of esculin, hydrolysis of gelatin and arginine dihydrazase. This species utilizes L-histidine and DL-malic acid as sole carbon and energy sources. Some strains utilize acetic acid, L-alanine, 3-hydroxy-butyrate, lactic acid, malonic acid and suberic acid while none of these strains utilize p-hydroxy-benzoate, citric acid, gluconic acid, propionic acid, L-serine or n-valeric acid. 18:1(w9c) is the major fatty acid and Q8 is the major quinone. DNA G+C content was 43-44 mol% as determined by HPLC. NIBH P2K6 isolated from seawater collected from the Japan Trench in the offing of Hachijo island,

Japan at a depth of 6,000m at 4°C , is used as a reference strain. This strain was deposited at the Institute for Fermentation, Osaka (IFO 16270). *Psychrobacter pacificus* strain NIBH P2K6 was deposited (IFO 16270) at the Fermentation Research Institute, 1-1-3, Higashi, Tsukuba-shi, Ibaraki-ken, Japan) dated May 21, 1999 (Accession No.: FERM BP-7106).

[0017]

Psychrobacter pacificus is a newly-discovered species and a number of its analogs have been isolated from seawater collected near the South Pole (Bowman et al., (1996) Int. J. Syst. Bacteriol. 46:841-848, Bowman et al., (1997) System. Appl. Microbiol. 20:209-215). Therefore, it will be appreciated that it is useful as an indicator organism with regard to the global circulation of deep-layer sea water. As one characteristic aspect of the global circulation of deep-layer sea water, it is known that deep-layer sea water in the Pacific Ocean flows steadily from the South Pole to the Japan Trench (Stommel, H.(1958) Deep-Sea Res. 5:80-82).

[0018]

16S rDNA having the base sequence represented by SEQ ID NO:1 can be obtained from *Psychrobacter pacificus* NIBH P2K6. Particularly, genomic DNA may be extracted from *Psychrobacter pacificus* NIBH P2K6 cells by any standard method, and 16S rDNA may be then amplified by PCR using appropriate primers (Lane, D.L. (1991) 16S/23S rRNA sequencing. In Nucleic Acid Techniques in Bacterial Systematics, pp.115-175. Edited by E. Stackebrandt, M.Goodfellow. West Sussex: John Wiley & Sons). An excess amount of primer and dNTP may be removed from the resulting PCR products and then the purified PCR products can be sequenced directly by cycle sequencing process using appropriate primers. The determined sequence is shown in SEQ ID NO:1.

[0019]

Based on the base sequence information of 16S rRNA gene from

Psychrobacter pacificus strain NIBH P2K6, a certain region in the base sequence specific to the bacterium may be extracted and a DNA probe may be prepared which enables molecular or cell level detection of the bacterium. Regions specific to the bacterium include one comprising nucleotide Nos. 458 to 476 of the nucleotide sequence of SEQ ID NO:1 (or nucleotide Nos. 469-487 in the corresponding *E. coli*. sequence. A DNA probe of 10-50bp in size, preferably 15-25bp in size, which corresponds to this region may be prepared. One preferable example may be a probe having the following base sequence:

• 5' TAATGTCATCGTCCCCGGG 3' (SEQ ID NO:2)

[0020]

A probe can be synthesized by phosphoramidite process (Beacage and Carruthers, Tetrahedron Lett. 22:1859-1862 (1981)) or triester method (Matteucci et al., J. Am. Chem. Soc. 103:3185 (1981)). Alternatively, a probe may be synthesized in an automatic synthesizer.

Further, a probe may be labeled with an isotope, a fluorescent label, DIG (digoxigenin) or the like. Examples of a label may include fluorescent colorants such as Cy5 (indodicarbocyanine), TRITC (tetramethyl rhodamine isothiocyanate) or FITC (fluorescein isothiocyanate), and haptens such as DIG (digoxigenin).

[0021]

A variety of processes for hybridization (e.g., Southern Blotting, Northern Blotting, Colony Hybridization or in-situ hybridization such as FISH; Fluorescence In Situ Hybridization) using a probe according to the present invention can be used to detect or identify *Psychrobacter pacificus*, *Psychrobacter glacincola* and analogs thereof. Analogs of *Psychrobacter pacificus* and *Psychrobacter glacincola* may include those contained in a database which do not have well-grounded identifications, *Psychrobacter glacincola* (AFO25555, PGU85879, PGU85878, PGU85877, PGU85876),

Psychrobacter immobilis (PIU85880), *Psychrobacter* sp. (PSU85874) and the like.

[0022]

One example of a process for detecting or identifying *Psychrobacter pacificensis* by using a probe having the nucleotide sequence of SEQ ID NO:2 will be described below.

A microorganism sample fixed by using, for example, paraformaldehyde is applied to a glass slide containing an organic film such as gelatin film formed thereon to immobilize the microorganism cells on the organic film. After dehydrating with ethanol or drying the microorganism cells at room temperature overnight, genomic DNA and RNA from the microorganism are allowed to hybridize to the DNA probe, and free- or incompletely bound-DNA probes are removed by washing. Microorganism cells are observed according to a conventional fluorescence microscopy procedure to detect fluorescence of the fluorescence-labeled DNA probe which has complementarily hybridized to nucleic acids from the subject microorganism. As a control, the same experiment may be performed using a non-specific DNA probe and a different microorganism belonging to another species. When a microorganism for which the DNA probe is targeted is used, it can be detected or identified since nucleic acids in cells will complementarily hybridize to the DNA probe thereby causing the cells to emit fluorescence.

The probe according to the present invention will not only enable species-specific detection of *Psychrobacter pacificus*, *Psychrobacter glacincola* and analogs thereof, but also will be very useful to enable a quicker and more accurate detection.

[0023]

[Example]

The present invention will be described in detail with reference to

the following examples which are included for purposes of illustration only and are not intended to limit the scope of the invention.

Example 1

Isolating *Psychrobacter pacificus* strains and sequencing 16S rRNA gene

Among 67 strains isolated from surface waters and the abyssal environment of the Japan Trench, 16 strains in total were selected. Eleven of the 16 strains were tentatively identified to be bacteria similar to *Moraxellaceae*. Different agar media such as 1/2 TZ based on an artificial seawater containing polypeptone and yeast extract (Maruyama, A et al., (1993) J. Oceanogr. 49, 353-367), Marine Agar (Difco; Detroit, MI, USA), Nitrient Agar (Difco) and the like were used to purify the strains. Each strain was incubated at 20°C and collected to obtain genomic DNA. A 1/2 TZ semi-solid agar medium containing 0.3% agar was used for storage at 4°C. Dry cells enclosed in a glass test tube was also stored at 4°C, for a long period of time.

[0024]

16S rRNA gene from those strains similar to *Moraxellaceae* isolated from the deep-sea water was determined by direct sequencing described below. Particularly, cells were collected by centrifugation, washed and resuspended in a TE buffer (10 mM Tris-HCl, 1 mM EDTA; pH 8.0). Cetyltrimethylammonium bromide (CTAB), phenol and chloroform/isoamyl alcohol (Murray, M.G. et al., (1980) Nucleic Acids Res 8, 4321-4325) were used to extract genomic DNAs by standard methods. In order to obtain almost-complete-16S rRNA gene, primers 27f and 1525r (Lane, D.L. (1991) 16S/23S rRNA sequencing. In Nucleic Acid Techniques in Bacterial Systematics, pp.115-175. Edited by E. Stackebrandt, M. Goodfellow. West Sussex: John Wiley & Sons) and the PCR cycle described in Maruyama, A. et al., (1997) Mar. Biol 128, 705-711 were used to perform PCR

amplification with Gene Amp PCR 9600 (Perkin-Elmer, Norwalk, Conn., USA). SUPREC-02 (Takara Shuzo Co., Ltd.) was used to remove an excess amount of primer and dNTP from the resultant PCR amplification products. An automatic DNA sequencer (ALFred; Pharmacia LKB, Sweden) was used to directly sequence the purified PCR products by cycle sequencing process according to the manufacturer's instruction using appropriate forward and reverse primers (Lane, D.L. (1991), *supra*). Particularly, the above-described primers were, *E. coli* numbering, 342r, 359f (5'-TCC TAC GGG AGG CAG CAG TG (SEQ ID NO:3); 20-mer), 519r, 803r (5'-CAT CGT TTA CGG CGT GGA C (SEQ ID NO:4); 19-mer), 821f (5'-GTC CAC GCC GTA AAC GAT G (SEQ ID NO:5); 19-mer), 1104r (5'-TTG CGC TCG TTG CGG GAC (SEQ ID NO:6); 18-mer), and 1111f (5'-GTC CCG CAA CGA GCG CAA (SEQ ID NO:7); 18-mer). Both strands of each fragment of 16S rDNA region were sequenced and ligated by using GENETYX software (version 8; Software Development Co., Ltd.). Except for the abyssal strains similar to *Moraxellaceae* of which 16S rDNA had been analyzed as mentioned above, 16S rDNA was extracted from other strains including the standard species *Moraxella lacunata* (ATCC 17967) as previously described in Maruyama, A. et al., (1997) (*supra*), amplified and subcloned. A multi-alignment program in CLUSTAL W (version 1.71; 44) was used to align these sequences. CLUSTAL W profile alignment option was used to align the sequences determined by the present inventors to known aligned sequences obtained from the rRNA www server ([http://rrna,uia,ac,be/](http://rrna.uia.ac.be/);45) of University of Antwerp. All locations which include a gap or gaps (i.e., undetermined or not-well-grounded sequences) were removed from the aligned data matrix. The nucleotide sequence of 16S rRNA gene from *Psychrobacter pacificus* NIBH P2K6 is shown in SEQ ID NO:1.

[0025]

Example 2

Result of Database Search for Probe Prepared

A probe sequence having the nucleotide sequence of SEQ ID NO:2 (named "Psypac469-487") was searched in RDP-DB (database) under conditions that permit up to 2 mismatches in the sequence. The sequences corresponding to Psypac469-487 were *Psychrobacter glacincola* (Pinhass et. Al)(0 mismatches), *Psychrobacter* sp., *Psychrobacter immobilis* and 4 strains involved in *Psychrobacter glacincola* , and *Psychrobacter* sp. Ant9 and *Moraxella* sp. Ant7 (Pinnhassi et. al.)(2 mismatches). These data indicate that Psypac469-487 can specifically detect strains similar to *Psychrobacter glacincola* and analogs thereof as well as *Psychrobacter pacificus*.

[0026]

These results indicate that the base sequence of the produced probe is complementary to the base sequence of *Psychrobacter pacificus*, *Psychrobacter glacincola* and analogs thereof, in a DNA database that encompasses all of the existing strains, thus showing that the probe may be very useful for species-specific detection of the two species.

[0027]

Example 3

Hybridization Test of Probe Prepared (1)

(1) Preparation of microorganism sample

Microorganisms (genus, species and strain) used herein are listed in Table 4. *Psychrobacter pacificensis* strains were isolated as microorganisms viable at 4°C from a deep-sea water sample collected from the Japan Trench in the offing of Hachijo island, Japan at a depth of 6,000m, but they could not be found at all in surface water (Maruyama et al. Marine Biology 128, 705-711, 1997). *Bacillus marinus* was cultured under an aerobic condition at 20°C using Marine Broth (Difco), and *Psychrobacter phenylpyruvicus* was cultured under an aerobic condition at

30°C using ATCC Culture Medium #4. The rest of the bacteria were cultured under aerobic conditions at 10-20°C using 1/2TZ liquid media (Maruyama et al., J. Oceanogr. 49, 353-367, 1993).

[0028]

Cultured microorganisms were fixed at 4°C overnight using a final concentration of 3% paraformaldehyde. The paraformaldehyde was first dissolved in 3×PBS (Phosphate Buffer Saline: 24g of NaCl, 0.6g of KCl, 0.72g of Na₂HPO₄, pH7.4) to a concentration of 15%, and then used in combination with samples (sample: paraformaldehyde=4:1).

[0029]

(2) Staining Sample with DNA Probe

The fixed microorganism sample was adsorbed on a Teflon-coated slide with a sample hole (diameter: 11mm) formed therein to which gelatin had previously been applied, and ethanol-dehydrated or dried at room temperature overnight. Next, 50 μl of a hybridization solution (0.9M NaCl, 50mM sodium phosphate buffer (pH7.0), 5mM EDTA, 0.5% SDS, Denhardt solution (final×1), 1.0mg/ml Poly(A)) was added to the sample hole. The glass slide prepared as described above was left to stand in a 50ml conical tube, loaded with a small amount of 3×PBS to prevent drying, at approximately 42°C for 30 minutes for prehybridization.

[0030]

Oligonucleotide DNA probes having the 5' end thereof fluorescence-labeled with Cy5, TRITC (tetramethyl rhodamine isothiocyanate) or FITC (fluorescein isothiocyanate) were prepared and each added to the above-described hybridization solution at an amount of 1ng probe/μl solution, and hybridization was performed at the temperature which is optimal for each oligonucleotide DNA probe for 4.5 hours. Unreacted oligonucleotide probe DNA was removed by washing the glass slide with wash solution (0.9M NaCl, 0.5mM sodium phosphate, 0.1% SDS,

pH=7.0) at the temperature which is suitable for each oligonucleotide DNA probe for 30 minutes.

[0031]

Oligonucleotide DNA probes used herein were Psypac469-487, a common probe among Domain Bacteria (previously Eubacteria) Eub338-355 (5'-GCTGCCTCCCGTAGGAGT (SEQ ID NO: 8), and a control probe, Cont (5'-GTGCCAGCAGCCGCGG (SEQ ID NO:9)).

[0032]

(3) Staining Sample DNA

After completion of hybridization, DAPI solution (final concentration of 5 μ g/ml) was added to the glass slide and left to stand at room temperature for 10 minutes to stain the DNAs present in the microorganism cells. After completion of the reaction, the glass slide was immersed in and washed with pure water for 15 minutes and then dried at room temperature.

[0033]

(4) Observation of Sample by Fluorescence Microscopy

An anti-color-degradation agent such as DABCO (diazabicyclooptane) solution (1 g/100 ml (10 ml PBS+90 ml Glycerol)) was added to the microorganism sample dried on the glass slide which was then covered with a cover glass and observed through a fluorescence microscope under oil immersion condition. Optionally, a fluorescent image for each colorant was captured by a cold CCD camera attached to the fluorescence microscope, and then analyzed. The results are shown in Table 4 below.

Table 4

Results of Utility Test of Probe by Fluorescence Microscopy According to
FISH

Strain	Fluorescence-labeled DNA probe				DAPI Staining
	Control	Pspac 469-487	Euba 338-355	Univ 1390-1407	
<i>Psychrobacter pacificensis</i> NIBH P2J2	x	○	○	○	○
<i>P. pacificensis</i> NIBH P2J3	x	○	○	○	○
<i>P. pacificensis</i> NIBH P2J13	x	○	○	○	○
<i>P. pacificensis</i> NIBH P2K6 (=IFO 16270)	x	○	○	○	○
<i>P. pacificensis</i> NIBH P2K18	x	○	○	○	○
<i>Psychrobacter glacincola</i> ACAM 483*	x	x	○	○	○
<i>Psychrobacter frigidicola</i> ACAM 304	x	x	○	○	○
<i>Psychrobacter immobilis</i> ATCC 43116	x	x	○	○	○
<i>Psychrobacter urativorans</i> ATCC 15174	x	x	○	○	○
<i>Psychrobacter phenylpyruvicus</i> ATCC 23333	x	x	○	○	○
<i>Pseudomonas aeruginosa</i> IFO 12689	x	x	○	○	○
<i>Vibrio parahaemolyticus</i> IFO 12711	x	x	○	○	○
<i>Bacillus marinus</i> ATCC 29841	x	x	○	○	○

* Previously described as *P. endoglaciecola* in the DNA database and later renamed for registration as a new species.

[0035]

These results show that the probes prepared herein can species-specifically bind to *Psychrobacter pacificus* and *Psychrobacter glacincola* but not to microorganisms belonging to any other genera, by *in situ* hybridization.

[0036]

[Advantage of the Invention]

Use of an oligonucleotide probe according to the present invention enables highly-sensitive and accurate, molecular- or cell-level detection of *Psychrobacter pacificus*, which is an indicator organism useful in understanding the circulation of deep-layer sea water. Since the use of a great number of microorganism samples is required for analysis of the behavior of deep-sea water and the evaluation of effects, enormous labor and time were required to perform complicated separating and culturing procedures as well as to classify and identify the microorganisms on land in conventional culture methods. However, use of the DNA probe according to the present invention, which comprises a base sequence, the species-specificity of which has been confirmed in an existing database, can provide quicker and easier detection or identification of *Psychrobacter pacificus*, *Psychrobacter glacincola* and analogs thereof.